

SK

ENVIROCHEM

Tableware Presoak

425 Whitehead Avenue
South River, NJ 08882
Phone—732-238-6700
Fax— 732-238-5590

GENERAL DESCRIPTION: SK TABLEWARE PRESOAK is a free-flowing liquid designed for presoaking silverware and stainless steel tableware. Its mild alkaline formula, including synthetic detergents and oxide removal agents, provides rapid penetration of food soils and discoloring films. SK TABLEWARE PRESOAK is completely safe on stainless steel and plastic.

PROPERTIES:

Chemical composition Synthetic soil cutting
and soil suspending agent; penetrating agents;
water conditioning agents.
Appearance Clear blue liquid
pH 11
pH @ 1% solution 9.4
Specific Gravity 1.15 (9.6 pounds per gallon)
Stability 1 Year @ Ambient Temperature
Wetting Ability Excellent
Foaming Moderate
Solubility Complete and Fast

INSTRUCTIONS FOR USE:

Use 2 oz. to 4 oz. of LIQUID SHINE per 1 gallon of water.
As a presoak for film removal on stainless steel tableware,
use same concentration and temperature.

RECOMMENDED RANGE OF USE CONCENTRATIONS:

Titrate with Institutional Alkalinity Test Kit. Use range is 0.1% to 0.3%, or 3 to 9 drops (each drop = .033%) using a 10 ml sample and 3 drops of Phenolphthalein Indicator. Select actual use concentration to handle water hardness and soil load.

NOTE: To properly titrate, adjust for raw water alkalinity by titrating H₂O just as you would a detergent solution. Subtract that number from total drop titration to derive intended detergent concentration.

Example: Raw water titrates 2 drops, total drop titration equals 8 drops. (8 - 2 = 6 x .033% = 0.2%)

HANDLING PRECAUTIONS:

**CAUTION: WARNING
KEEP OUT OF REACH OF CHILDREN.**

Highly alkaline liquid. Contact with skin and eyes may cause severe irritation. Harmful if swallowed. Flood with water if splashed in eyes. Flush with water for 15 minutes and consult physician immediately. If ingested, drink large amounts of milk, milk of magnesia or gelatin, or if these are not available, drink large amounts of water. DO NOT induce vomiting. Obtain prompt medical attention.

PACKAGING:

5 gal plastic pail - 54 lbs. net / 57 lbs. gross
15 gal plastic drum - 162 lbs. net / 168 lbs. gross
55 gal plastic drum-583 lbs. net/600 lbs. gross

Seller makes no warranty, expressed or implied, concerning the use of this product other than indicated on the label. Buyer assumes all risk of use and/or handling of this material when such use and/or handling is contrary to label instructions.